

High Resolution Elemental Mapping of Biological Samples by High Sensitivity, Fast-Scanning LA-ICP-MS

Lothar Rottmann¹, Julian D. Wills¹, Bo Yan², Catherine E. Costello², Lee E. Goldstein², Dhinesh Asogan³
¹Thermo Fisher Scientific (Bremen) GmbH, 28199 Bremen, Germany; ²Boston University School of Medicine, Boston, MA 02218, USA; ³Teledyne CETAC Technologies, Omaha, NE 68144, USA

Overview

Purpose: To demonstrate the use of LA-ICP-MS for high spatial resolution metallomic imaging mass spectrometric analysis of biospecimens relevant to medical research on traumatic brain injury.

Methods: Samples were prepared and subsequently analyzed by LA-ICP-MS. Metallomic distribution images were constructed using in-house built MATLAB scripts.

Results: Elemental and isotopic images for a series of biological samples were produced and evaluated with results from optical imaging.

Introduction

We have developed a technique – metallomic imaging mass spectrometry (MIMS) tissue mapping – to localize and quantify neuroinflammation and disruption to the blood-brain barrier (BBB) following experimental traumatic brain injury (TBI). This novel method provides a wide detection range and permits single-run, multi-elemental and/or isotopic mapping in tissue specimens with variable elemental/isotopic concentration and composition .

Methods

Instrumentation

Mouse brain samples were sampled using a Teledyne CETAC LSX-213 G2+™ laser ablation system equipped with the HelEx II active 2 volume laser cell for optimum washout. Laser parameters:

- Spot size: 20 µm
- Shot frequency: 20 Hz
- 100% laser energy
- Single line scan
- Pause between lines: 30 s.

The ablated material was carried in a helium stream to either a quadrupole (Thermo Scientific™ iCAP™ Q) or high resolution (Thermo Scientific™ ELEMENT XR™) ICP-MS for analysis.

FIGURE 1. Combined LSX-213 G2+ / iCAP Q LA-ICP-MS system



Results

Mapping the metallome of the mouse hippocampus by MIMS (ELEMENT XR ICP-MS).

FIGURE 2. Mouse brain, Nissl stain. Red box, hippocampus and dentate gyrus, critical brain regions for memory formation. Allen Institute for Brain Science.




FIGURE 3. Mouse hippocampus and dentate gyrus showing areas of zinc enrichment by Timms stain. (Frederickson CJ, et al., *J Nutr* 130:1471, 2000).



FIGURE 4. ^{56}Fe MIMS map (20 μm resolution)




FIGURE 5. ^{65}Cu MIMS map (20 μm resolution)




FIGURE 6. ^{70}Zn MIMS map (20 μm resolution)




FIGURE 7. Composite ^{56}Fe , ^{65}Cu , ^{70}Zn MIMS map with elemental concentrations at a selected point (CA1 subfield), boxed region (dentate gyrus), and segmentation area (CA3 subfield).





FIGURE 8. $^{157}\text{Gd-MIMS}$ (iCAP Q ICP-MS) mapping of the mouse brain after experimental traumatic brain injury (TBI).



Note the correspondence between the MRI-detectable enhancement in the left perirhinal cortex (arrows, top two rows) and Gd accumulation in the $^{157}\text{Gd-MIMS}$ map (arrow, bottom row). The $^{157}\text{Gd-MIMS}$ map provides direct quantitative evidence of anatomically localized BBB disruption following impact TBI. The closed-head concussive impact injury induced only minor contusion (third row). These MIMS results provide clinical insights into the pathophysiological mechanisms of concussion and acute neurotrauma (taken from Tagge et al., in preparation).

Conclusion

- MIMS analysis enables high-throughput, high spatial resolution, multi-elemental analysis of biological samples with definitive elemental/isotopic identification.
- Metallomic distribution provides clinical insights in the blood brain barrier dysfunction and neuroinflammation that are not provided by traditional techniques.

Acknowledgement

Research support: DoD TATRC Award W81XWH-13-1-0263 (LEG), NIH P41 GM104603 (BY & CEC), BUADC Pilot Grant 9500-304-584 (BY).

www.thermofisher.com

© 2016 Thermo Fisher Scientific Inc. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. LSX-213 G2+ is a trademark of Teledyne Instruments Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. For research use only, not for diagnostic purposes. P043316 First presented at the Winter Conference on Plasma Spectrochemistry, Tucson, AZ, January 2016

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 10 3292 200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
 SCIENTIFIC

A Thermo Fisher Scientific Brand